

4D Hyaluronic Acid (Hymagic™-4D)

INCI name: Sodium Hyaluronate crosspolymer, Sodium Hyaluronate, Sodium Acetylated Hyaluronate, Hydrolyzed Sodium Hyaluronate, Pentylene Glycol, Ethylhexylglycerin, Aqua

Specification

1. Characters	A clear transparent, viscous solution
2. Colour	Colorless to slight yellow
3. Odour	Slight characteristic odour
4. Identification	
Reaction of uronic acid	Positive
5. pH	5.0~7.0
6. Transparency	≥97.0%
7. Kinetic viscosity	≥500 mPa•S
8. Heavy metals	≤10ppm
9. As	≤2ppm
10. Bacteria counts	≤100cfu/g
11. Mold & Yeast	≤10cfu/g
12. <i>Pseudomonas aeruginosa</i>	<i>Negative</i>
13. <i>Staphylococcus aureus</i>	<i>Negative</i>
14. <i>Escherichia coli</i>	<i>Negative</i>

【Storage】 Kept airtight in room temperature, protected from light.

【Package】 1kg/bottle, 20kg/barrel

【Validity】 Two years.

Recommend usage: 1%-10%